FLACONE 10 ml - 10 ml MULTIDOSE BOTTLE

Composizione/Composition: Purified Water, *Sea Water (28-30%), Jacobaea Maritima Extract (0,01%), Calendula Officialis Flower Water (0,10%), Centaurea Cyanus Flower Water (0,10%), Helichrysum Italicum Flower Water (0,10%), Euphrasia Officialis Flower/Leaf/Stem Water (0,10%), Hydroxypropylmethylcellulose (0,20%), Sodium Hyaluronate (0,20%). *Isotonic

IT ISOMAR® Occhi Gocce Oculari

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EN ISOMAR® Occhi Eye Drops

INDICATIONS FOR USE: **Isomar Occhi** provides relief from redness, irritation, burning, itching and the sensation of foreign bodies, as well as from tearing and fatigue induced by factors such as: exposure to wind and sunlight, dust, smog, smoke, dry heat, air-conditioned environments, prolonged exposure to computer screens and contact lens wear. **Isomar Occhi** contributes significantly in alleviating the discomfort caused by pollen coming into contact with the conjunctiva by relieving allergic symptoms. **Isomar Occhi** helps with dry eyes and can be used with all types of contact lenses. Application of **Isomar Occhi** prior to fitting contact lenses is recommended for an enhanced lubricating and soothing action.

INSTRUCTIONS FOR USE: to use the product correctly, follow the instructions below:

- open the package, remove the protective seal and cap when administering the first dose;
 - tilt the head and hold the vial upside down over the eye;
 - gently hold the lower eyelid down with one finger and gradually apply more pressure with the fingers on the bottle until a drop appears;
 - apply 1-2 drops into the conjunctival sac every 4-5 hours;
 - close the eye for a few seconds;
 - close the bottle immediately after application. The solution can be reused up to 6 months after the first opening. If required, the procedure can be repeated for the other eye.
- The cap used on **Isomar Occhi** vials is specifically designed to guarantee the sterility of the product during use and for optimal dispensing. It is recommended that a closed vial not be pressed directly to prevent any drops from escaping from the grille of the cap. As the solution in the bottle empties, it may be necessary to apply more force.

INTENDED PATIENT POPULATION: the medical device is intended for users with no specialist knowledge by observing the manufacturer's instructions for use. The medical device is intended for use in home environments. The user is provided with all the necessary information for the correct and safe use of the device.

CONTRAINDICATIONS: the medical device is not indicated for individuals with known hypersensitivity to the ingredients in the formulation.

WARNINGS: carefully read the information provided by the manufacturer before use. Do not exceed the recommended dose. Discontinue treatment if any undesirable effect occurs. Do not use in conjunction with other topical treatments. Do not swallow. The product is intended for ophthalmic use only. Store inside the box in a dry location, away from sources of light and heat, at a temperature between 5°C and 25°C. Do not pierce or tamper with the packaging. Do not use the product if the package is opened or damaged. Do not use the product after the expiry date. The expiry date refers to the unopened and properly stored product. The product is intended for adult patients. There is no specific data available regarding the use of the product on children and pregnant women. Keep out of sight and reach of children. Wash hands thoroughly before applying the solution and do not touch the tip of the dispenser with the fingers. When administering the solution, avoid touching the conjunctival membrane with the tip of the vial. The solution must be used up to 6 months after the first opening of the vial. For correct storage, accurately close the bottle after use. Do not dispose of the vial in the environment after use. Any severe accident resulting from the use of the medical device must be immediately reported to the manufacturer and the competent authority of the member state in which the event occurred.

DOSAGE: 1-2 drops per eye every 4-5 hours. We recommend alternating treatment cycles with suspensions of 7 days. Each cycle should not exceed 30 days of treatment.

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